

Tobii Dynavox Compass

Text Messaging Quick Start Guide

Gather your Devices

In order to use text messaging with Tobii Dynavox Compass, you will need an Android-based Tobii Dynavox T-series device (T7, T10, T15) and an Android smart phone.



Your Tobii Dynavox device and the Android smart phone will communicate via Bluetooth, so the phone must remain on and within approximately 32 feet (10 meters) of your Tobii Dynavox device for text messaging to work from Compass.

Enable Bluetooth

Follow the steps below on both your devices

1. Select **Settings** to open the Android settings menu.
2. Set Bluetooth to the **ON** position.

Download and Install the Tablet Talk Software



If you purchased your Tobii Dynavox device after April of 2016, you will already have the updated Tablet Talk software installed on your Tobii Dynavox device. Skip to *Install Tablet Talk on your Android Phone* to download and install the necessary software to your phone.

Update Tablet Talk on your Tobii Dynavox Device

1. On your Tobii Dynavox device, open the web browser and navigate to http://www.mytobiidynavox.com/Content/downloads/TT_180_TABLET_X2_836.apk
2. Save the file to your device.
3. Launch the file and follow the installer instructions to update Tablet Talk.

Install Tablet Talk on your Android Phone

1. On your phone, open the web browser and navigate to http://www.mytobiidynavox.com/Content/downloadsTT_180_X_836.apk
2. Save the file to your phone.
3. Launch the file and follow the installer instructions to install Tablet Talk.

Pair Your Devices

Once you have enabled Bluetooth and installed Tablet Talk on both of your devices, you can then pair your phone with your Tobii Dynavox device via Bluetooth.

Support for Your Tobii Dynavox Device


Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.


Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

1. *On your device:* The Bluetooth Pairing Wizard should have opened after updating Tablet Talk. If the Bluetooth Pairing Wizard is not open, select the **wifi/Bluetooth** icon in the lower left corner of Tablet Talk, then choose **Setup**.



2. Select **Next**.
3. *On your phone:* Select **Enable Discovery**.



Support for Your Tobii Dynavox Device


Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.


Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

4. Select **Allow** to make your phone visible for Bluetooth discovery.



5. *On your device:* Select **Next**.
6. Find and select **your phone's name** in the list. A connect prompt will open.



Support for Your Tobii Dynavox Device


Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.


Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

7. Select **Okay**.



8. *On both your phone and your device:* You will see a passkey prompt on both your Tobii Dynavox device and your phone. Verify that the passkey matches and, if so, select **Pair** on both devices.



9. Your devices should now be paired and the Bluetooth icon in the Tablet Talk app will be green.



Send and Receive Text Messages

Send Text Message functionality has been added to the following English pagesets: All Access, Core First, Masterpage, Nav-Bar, Navigator, and Stroke and Brain Injury. Devices shipped after April 2016 will have these updated pagesets. All other users can download the updated pagesets from myTobiiDynavox.com.

Send a Text Message

1. *On your Tobii Dynavox device:* Launch **Compass**.
2. Select **Keyboard**.
3. In the message window, enter the text that you would like to send.



Compass text messaging only supports sending simple text messages. It cannot send symbols, pictures, or sounds via text message.

Support for Your Tobii Dynavox Device


Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.


Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact


4. Select **Send**.



5. Select **Apps**.



6. Select **Send Text Message**.



Support for Your Tobii Dynamox Device


Get Help Online

See the product-specific Support page for your Tobii Dynamox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynamox.com or www.myTobiiDynamox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynamox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynamox.com/contact


- If you have already entered the desired contact, skip to step 12. If you don't see the desired contact in the list, select **Create New Contact**.



- Type the contact name, then select **Enter Name**.



- Select **mobile phone number**, enter the phone number for the contact, then select **Done**.



- If desired, choose a profile picture for the contact.

Support for Your Tobii Dynavox Device


Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.


Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

11. Select **Save**.




12. Select the desired contact to send the text to that person.



Receive a Text Message

1. When you receive a text message, it will appear as a popup on top of the Compass app.



Contacts created in Compass are not replicated in the Tablet Talk app automatically. To see contact names and photos on incoming text messages, the contact information must be manually entered in the Tablet Talk app.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Reply to a Text Message

- a. People using any access method other than scanning can reply to the text directly in the Tablet Talk popup.
- b. People using the scanning access method must select the **app switcher** in the upper right corner, choose **Compass**, then compose a new text in Compass if they wish to reply.



Text messages received while the Tobii Dynavox device is off will not appear as a popup. They can, however, be viewed in the Tablet Talk app.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact