

Release Notes **Tobii Dynavox Windows Control**

Introduction

Tobii Dynavox Windows Control is a software solution for controlling Windows applications using gaze with Tobii Dynavox eye trackers.

System Requirements

- Windows 7, 8.1, or 10
- .NET Framework 4.5

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 2.0.2

Release date: 2017-10-16

New Features

- Added basic statistics recording (TGP-594)
- Added Korean keyboard and word prediction (TGP-1524)
- Added Sticky Move Cursor, Sticky Drag & Drop, Sticky Right Click and Sticky Double Click tasks (TGP-1981)
- Added a video option to use as a stimuli for Easy calibration (TGP-2360)
- Added video calibration stimuli samples (TGP-2483)
- Added a notification to inform users when Windows Eye Control is running at the same time with Windows Control (TGP-2484)

Improvements

- Made it possible to cancel out Display Setup and Calibration notifications (TGP-1467)
- Updated the Track Status visual design (TGP-1560)
- Fixed a bug where eye trackable regions had memory leak problems (TGP-1572)
- Made it possible to switch back to the previously selected calibration profile for the first 15 seconds while a new one is selected (TGP-1897)
- Fixed a bug where in some cases Windows Control keyboard was shown behind the Windows Control taskbar (TGP-1977)
- Made notification windows visible on top of Windows Control Settings (TGP-2135)
- Fixed a bug where a tertiary selection region remained on the previous screen after changing the eye trackable display (TGP-2141)
- Fixed a bug where the Windows Control taskbar didn't change a side when Charms were shown by Shortcuts in Windows 8 (TGP-2159)
- Fixed a bug where eye trackable regions were inconsistently cropped when they were covered by windows (TGP-2199)
- Fixed a bug where the eyes in the Windows Control taskbar didn't disappear when they weren't detected by the eye tracker (TGP-2223)
- Fixed a bug where showing tertiary selection regions using mouse would ignore "Show Additional Tasks" setting (TGP-2233)
- Fixed a bug where the eyes in the Windows Control taskbar and in the Track Status didn't disappear when the eye tracker was disconnected (TGP-2235)
- Changed the time limit for diagnostics information to be collected only for the last month period (TGP-2397)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

- Made it possible to set the bracket configuration in Display Setup notification window for relevant eye trackers (TGP-2449)
- Fixed a bug where the mounting mode wasn't updated in Windows Control Settings when the eye tracker was disconnected (TGP-2450)
- Updated EyeMobile Plus configuration measurements for the predefined screen planes (TGP-2451)
- Updated the tutorial with additional information about Sticky tasks (TGP-1981)
- Fixed a bug where continuous activations of eye trackable region didn't re-evaluate its visible bounds (TGP-2497)
- Improved the Windows Control performance by reducing the amount of WMI queries (TGP-2455)
- Fixed a bug where Tobii Eye Tracking Core continuously showed its context menu while Eye Assist Engine was running at the same time (TGP-2501)
- Fixed a bug where the scroll marker was shown on a wrong screen in cases when the primary monitor wasn't set as eye trackable (TGP-2504)
- Made it possible to step through Easy calibration using mouse clicks or by tapping on the screen on touch devices (TGP-2508)
- Changed the maximum timeout to 10 seconds for clients to register to Eye Assist Engine (TGP-2515)
- Fixed a bug where Windows Control Settings were shown in-between calibration and test calibration windows (TGP-2516)
- Added localization of the new text (TGP-2519)
- Made it possible to close the Windows Control keyboard using Enter key when it is shown in Name Profile window in calibration process (TGP-2522)
- Fixed a bug where eye trackable regions were shown behind notification windows (TGP-2524)
- Fixed a bug where clicking Esc keyboard key in Name Profile window in calibration process would close the calibration but not the Windows Control keyboard (TGP-2539)
- Updated the way of how the Windows Control keyboard state is detected (TGP-2540)
- Fixed a bug where calibration process crashed when the eye tracker was disconnected (TGP-2546)
- Fixed a bug where Virtual Remote Drivers application didn't detect its previous installation (TGP-2562)
- Fixed a bug where changing "Limit Eye Gaze buttons only to Tobii Dynavox applications" setting would disable eye trackable region activations (TGP-2573)
- Fixed a bug where tertiary selection region wasn't positioned properly (TGP-2616)
- Fixed a bug where Windows Control crashed after docking the computer in sleep mode and then waking it up (TGP-2525)
- Updated the bundled version of Virtual Remote to 1.0.4 (TGP-2608)
- Updated the bundled version of Update Notifier to 1.5.3 (TGP-2557)
- Updated the bundled version of Tobii Eye Tracking Core to 2.10.20 (TGP-2608)

Release Version 2.0.1

Release date: 2017-07-10

New Features

- Support for EyeMobile Plus (TGP-2084, TGP-2373)

Improvements

- Installers installed by the Windows Control bundle can now be repaired through Windows (TGP-1666)
- Made it possible to switch to sticky left click using mouse (TGP-1929)
- Corrected Windows Control Configuration Guide welcome page button's font (TGP-2032)
- Fixed inconsistent highlighting in the manage profile view in Settings (TGP-2099)
- Fixed a bug where the incorrect calibration type was shown if the serial number of the eye tracker could not be retrieved (TGP-2113, TGP-2210)
- Improved look and feel of the track status window in it's smallest size (TGP-2146)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

- Fixed a bug where the Shortcuts window would close when pausing Windows Control from the Controller window (TGP-2166)
- Fixed an incorrect Portuguese translation for the eye tracking disabled notification (TGP-2183)
- Added pause and resume to the tray icon's context menu (TGP-2193)
- Fixed a bug where the "More Info" button in Settings was not possible to click using mouse when gaze interaction was paused (TGP-2196)
- Fixed a bug where the feedback dot would show in the Windows Control zoom window when using secondary sticky left click (TGP-2200)
- Fixed an error in the Windows Control shutdown process (TGP-2201)
- Fixed a bug where the Controller window could end up outside of all screens when changing eye tracked screen (TGP-2204)
- Fixed a typo, "Windows", for Norwegian (TGP-2205)
- Fixed a bug with the distance meter of the track status when using PCEye Mini (TGP-2221)
- Updated the bundled version of Tobii Eye Tracking to 2.8.2 (TGP-2230)
- Renamed the 'Next' button to 'Start' in the Calibration Name Profile (TGP-2240)
- Fixed a bug where regions with negative dwell times could crash Region Interaction (TGP-2244)
- Fixed a bug where the Windows Control taskbar could become mispositioned if another appbar is docked to top (TGP-2249)
- Updated font on the Controller window's button (TGP-2264)
- Fixed incorrect Chinese translation of "Activation Method" (TGP-2265)
- Fixed a bug where the Windows Control keyboard would get stuck on the screen when locking and unlocking Windows while the keyboard was open (TGP-2290)
- Fixed a bug where the system information file would not be created if any of its intended content could not be accessed (TGP-2307)
- Fixed cut off text in the notification for Tobii Service not running (TGP-2309)
- Fixed a deadlock in Region Interaction when opening and closing the Windows Control keyboard from other applications (TGP-2310)
- Fixed a bug where there could be a lingering feedback dot after turning Windows Control off (TGP-2316)
- Changed the default dwell time to 1100 milliseconds (TGP-2331)
- Updated the taskbar task selection UI in Settings to better fit translated texts (TGP-2337)
- Updated the tutorial UI to better fit translated texts (TGP-2350)
- Made the Test Eye Tracking view topmost (TGP-2355)
- Fixed inconsistent capitalization of text on the tiles of the Setting application's start page (TGP-2357)
- Fixed a bug in the German Configuration Guide where "Konfigurationsassistent" would split to two lines (TGP-2393)
- Updated the bundled version of Virtual Remote to 1.0.1 (TGP-2417)
- Updated the bundled version of Update Notifier to 1.5.1 (TGP-2444)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact