

Release Notes – TET Server 2.9.16

Compatibility issues

- TET Server 2.7.13 and later is *not* backwards compatible with older versions. Falling back to version 2.6.3 or older needs special instructions from Tobii Technology.

Tobii Eye Tracking Server 2.9.16

New or changed functionality

- **Automatic detection of product**
The TET Server automatically detects which product is connected. The product configuration is stored per unit instead of per product family.
- **Infant AddOn will NOT be marked as connected** in the tray icon after a TET Server upgrade. For customers using the Infant AddOn please:
 1. Verify that the eye tracker is connected and detected in the tray menu.
 2. Switch Infant AddOn to *connected* from the tray menu.Attention!! **This setting is now individual** for a specific eye tracker. The computer will remember what eye trackers it has been used with and it will remember whether or not an Infant AddOn was connected last time this eye tracker was used.
- **The x50 configuration will be set to default.** Those who have an x50 setup and wish to upgrade should:
 1. Save the current x50 configuration to a file. This can be done from the x50 config tool.
 2. Upgrade to TET Server 2.9.16.
 3. Connect the Tobii x50 and verify from the Tobii tray meny that it is detected.
 4. Open the x50 config tool and load the old configuration.Attention!! **This setting is now individual** for a specific eye tracker. Each Tobii x50 unit has to be configured separately.
- **Old calibrations cannot be improved**
Calibrations done using a TET Server older than 2.9.9 cannot be improved using TetServer 2.9.9 or newer. The calibrations can still be loaded.
- **Automatic improvement of the calibration**
When calibrating, the TET Server will automatically search for bad data and remove these. That might make the calibration process a little slower.