

## Release Notes **Gaze Point**

### Introduction

The Tobii Dynavox Gaze Point Software is a mouse move and dwell click application to be used with third party software.

### System Requirements

- Windows 7, 8, 8.1 or 10
- .NET Framework 4 Full

**Customer support:** For more information or support, please visit [www.TobiiDynavox.com](http://www.TobiiDynavox.com).

### Release Version 1.1.1

**Release date:** 2016-02-22

#### Improvements

- Fixed a bug where the mouse cursor would move very slowly when running on Windows 7 touch devices (TGP-391)
- Corrected Support and FAQ links (TGP-481, TGP-508)
- Fixed a bug where the positioning guide prior to calibration would get stuck and only 'Cancel' button would work to continue the process (TGP-679)
- PCEye Update Notifier 1.4.2 included (TGP-825)

### Release Version 1.1.0

**Release date:** 2015-09-15

#### New Features

- Support for Windows 10 (TGP-412)
- 2-point calibration with all supported eye trackers, and also 1-point calibration with PCEye Explore (TGP-251)
- Gaze Point is localized and support two new languages, German and French (TGP-418)
- Gaze Point can be updated through PCEye Update Notifier (TGP-422)

#### Improvements

- Support Tobii EyeX Engine 1.4.0 version (TGP-419)
- Upgrade to the latest eye tracker firmware (TGP-424)
- Replaced some calibration videos (TGP-426)
- Improved the position guide in the calibration flow (TGP-430)
- Fixed a bug where Gaze Point had the wrong link to the support website (TGP-414)

### Release Version 1.0.1

**Release date:** 2015-03-10

### Support for Your Tobii Dynavox Device

#### Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: [www.TobiiDynavox.com](http://www.TobiiDynavox.com) or [www.myTobiiDynavox.com](http://www.myTobiiDynavox.com).

#### Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit [www.TobiiDynavox.com/contact](http://www.TobiiDynavox.com/contact)



## Improvements

- Explicitly block non-supported Eye Trackers. (TGP-385)

## Support for Your Tobii Dynavox Device

### Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: [www.TobiiDynavox.com](http://www.TobiiDynavox.com) or [www.myTobiiDynavox.com](http://www.myTobiiDynavox.com).

### Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit [www.TobiiDynavox.com/contact](http://www.TobiiDynavox.com/contact)