

Release Notes **Gaze Point**

Introduction

The Tobii Dynavox Gaze Point Software is a mouse move and dwell click application to be used with third party software.

System Requirements

- Windows 7, 8.1, or 10
- .NET Framework 4.5

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 2.0.2

Release date: 2017-10-16

New Features

- Added basic statistics recording (TGP-427)
- Added a video option to use as a stimuli for Easy calibration (TGP-2360)
- Added video calibration stimuli samples (TGP-2483)
- Added a notification to inform users when Windows Eye Control is running at the same time with Gaze Point (TGP-2484)

Improvements

- Made it possible to cancel out Display Setup and Calibration notifications (TGP-1467)
- Updated the Track Status visual design (TGP-1560)
- Fixed a bug where eye trackable regions had memory leak problems (TGP-1572)
- Made it possible to switch back to the previously selected calibration profile for the first 15 seconds while a new one is selected (TGP-1897)
- Made notification windows visible on top of Gaze Point Settings (TGP-2135)
- Fixed a bug where eye trackable regions were inconsistently cropped when they were covered by windows (TGP-2199)
- Fixed a bug where the eyes in the Track Status didn't disappear when the eye tracker was disconnected (TGP-2235)
- Changed the time limit for diagnostics information to be collected only for the last month period (TGP-2397)
- Fixed a bug where object leakage occurred while closing Gaze Point using gaze (TGP-2428)
- Made it possible to set the bracket configuration in Display Setup notification window for relevant eye trackers (TGP-2449)
- Fixed a bug where the mounting mode wasn't updated in Gaze Point Settings when the eye tracker was disconnected (TGP-2450)
- Updated EyeMobile Plus configuration measurements for the predefined screen planes (TGP-2451)
- Fixed a bug where Gaze Point didn't restart after Gaze Point Settings were closed (TGP-2465)
- Fixed a bug where continuous activations of eye trackable region didn't re-evaluate its visible bounds (TGP-2497)
- Made it possible to step through Easy calibration using mouse clicks or by tapping on the screen on touch devices (TGP-2508)
- Changed the maximum timeout to 10 seconds for clients to register to Eye Assist Engine (TGP-2515)
- Fixed a bug where Gaze Point Settings were shown in-between calibration and test calibration windows (TGP-2516)
- Added localization of the new text (TGP-2519)
- Fixed a bug where eye trackable regions were shown behind notification windows (TGP-2524)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

- Fixed a bug where calibration process crashed when the eye tracker was disconnected (TGP-2546)
- Updated the bundled version of Update Notifier to 1.5.3 (TGP-2557)
- Updated the bundled version of Tobii Eye Tracking Core to 2.10.20 (TGP-2608)

Release Version 2.0.1

Release date: 2017-07-10

New Features

- Support using switch input to click (TGP-866)
- Support PCEye Mini, PCEye Plus and EyeMobile Plus (TGP-960)
- New Gaze Point Setting application (TGP-995)
- New Gaze Point Configuration Guide (TGP-1636)

Improvements

- Gaze Point works with region interaction (TGP-532, TGP-587)
- Updated user interface of the main controller window (TGP-1779, TGP-1198)

Release Version 1.1.2

Release date: 2016-08-15

Improvements

- Fixed problem with repeated "Eye Tracker connected" messages when using PCEye Mini through updated EyeX version (TGP-1032)
- PCEye Update Notifier 1.4.3 included (TGP-1096)

Release Version 1.1.1

Release date: 2016-02-22

Improvements

- Fixed a bug where the mouse cursor would move very slowly when running on Windows 7 touch devices (TGP-391)
- Corrected Support and FAQ links (TGP-481, TGP-508)
- Fixed a bug where the positioning guide prior to calibration would get stuck and only 'Cancel' button would work to continue the process (TGP-679)
- PCEye Update Notifier 1.4.2 included (TGP-825)

Release Version 1.1.0

Release date: 2015-09-15

New Features

- Support for Windows 10 (TGP-412)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

- 2-point calibration with all supported eye trackers, and also 1-point calibration with PCEye Explore (TGP-251)
- Gaze Point is localized and support two new languages, German and French (TGP-418)
- Gaze Point can be updated through PCEye Update Notifier (TGP-422)

Improvements

- Support Tobii EyeX Engine 1.4.0 version (TGP-419)
- Upgrade to the latest eye tracker firmware (TGP-424)
- Replaced some calibration videos (TGP-426)
- Improved the position guide in the calibration flow (TGP-430)
- Fixed a bug where Gaze Point had the wrong link to the support website (TGP-414)

Release Version 1.0.1

Release date: 2015-03-10

Improvements

- Explicitly block non-supported Eye Trackers. (TGP-385)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact