

Release Notes **Tobii PCEye Go**

Introduction

The Tobii PCEye Go Software is a software solution for setting up and configuring use of Gaze Interaction with the Tobii PCEye Go eye trackers.

System Requirements

- Windows Vista (SP2) 7, 8, 8.1 or 10
- .NET Framework 4 Full
- DirectX 9.0c

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 2.3.0

Release date: 2015-09-15

New Features

- Support for Windows 10
- Version 2.7.0 of Gaze Interaction Software included
- Version 1.4.0 of Update Notifier included
- Installer for Gaze Selection word prediction language packs included

Improvements

- Name change from Tobii to Tobii Dynavox. (CX-623)
- Improved error handling when the PCEye Configuration Guide is running without Gaze Interaction Software installed. (CX-639)

Release Version 2.2.4

Release date: 2014-09-20

New Features

- Version 2.6.1 of Tobii Gaze Interaction Software included
- Version 1.3.1 of PCEye Update Notifier included
- Updated the user manuals with illustrations of new magnetic mount (CX-613)

Improvements

- Resolved issue where a non-executable file was set to run at Windows startup (CX- 610)

Release Version 2.2.3

Release date: 2014-05-16

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Improvements

- Updated the user manuals. (CX-601)
- Added device data for Microsoft Surface 2. (CX-590)
- Added device data for Dell Venue 11 Pro 5130. (CX-590)
- Added device data for Dell Venue 11i Pro 7130/7139. (CX-590, CX-611)

Release Version 2.2.2

Release date: 2013-12-02

Improvements

- Updated the user manuals. (CX-592)
- Added device data for Microsoft Surface Pro 2. (CX-593)

Release Version 2.2.1

Release date: 2013-11-06

Improvements

- Updated the pictures in the Installation Guide to match the new mounting plate. (CX-586)
- Corrected the metadata of the self extracting installer. (CX-576)

Release Version 2.2.0

Release date: 2013-10-11

New Features

- New version of Tobii Gaze Interaction Software; 2.5.0.
- New version of EyeTrackerBrowser; 3.0.83.
- Update Notifier. This application notifies the user of available software updates and installs these updates.

Improvements

- Improved values for the Positioning Guide track box. (CX-440)
- Fixed a bug where PCEye Configuration guide could launch Windows Control unexpectedly. (CX-538)

Release Version 2.1.0

Release date: 2013-09-06

New Features

- Installs Tobii Gaze Interaction Software 2.4.0.
- Support for EyeMobile.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Improvements

- .NET 3.5 now included in PCEye Installation Guide. (CX-399)
- First page next button in PCEye Installation Guide is now disabled until the user selects a product. (CX-528)
- Improved performance of track status in PCEye Configuration Guide. (NESDK-1082)
- The Configuration Guide now only shows one confirmation regardless of how many time the save screen size button is pressed. (PCEYE-31)
- The Configuration Guides Screen and Calibration pages are now enabled/disabled properly with attached/removed EyeTracker. (PCEYE-47)
- The Installation- and Configuration Guides are now always possible to exit when they are complete. (PCEYE-48)
- Tobii Gaze Interaction Settings might crash if the process is started as another Windows user than the owner of the active session. (CX-439, CX-464)
- Server communication failure if Tobii Gaze Interaction Software is uninstalled prior to PCEye software while Notifier.exe is running. (CX-509)
- The PCEye Installation Guide window is too large when using 125% scaling on devices with low screen resolution. (CX-531)

Release Version 2.0.1

Release date: 2013-07-08

New Features

- Installs Tobii Gaze Interaction Software 2.3.1.

Improvements

- Support for the old PCEye in the new PCEye Configuration Guide. (CX-435)
- The USB installation stick now only shows one installation option for the user. (CX 510, PCEYE-8)
- Manuals are accessible on the USB installation stick prior to installation. (PCEYE-9)
- Updated the text "remove USB now please" to "Please remove the Tobii PCEye USB flash drive.". (CX-511)

Release Version 1.3.0

Release date: 2012-08-29

Improvements

- The Installation Guide, Configuration Guide and Tobii Notifier are all upgraded to use .NET 4 (Internal id: CX-329)
- Upgrade to use Tobii Eye Control Suite 2 (Internal id: CX-328)
- The default Windows Control method for PCEye is now Tobii Gaze Selection (Internal id: CX-330)
- The texts in the PCEye Installation and Configuration Guides become invisible when using "High Contrast Black" color scheme (Internal id: CX-314)
- The input fields for centimeters and inches on the screen measurement page in the PCEye Configuration Guide interact with each other (Internal id: CX-332)
- It is possible to make the PCEye Configuration Guide accept the screen size 0 x 0 cm (Internal id: CX-334)
- Tobii Notifier does not start up automatically after PCEye Configuration Guide (Internal id: CX-338)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact