

Communication device
for increased control and independence



Tobii C12 Communication Device with Eye Control



Communicate, access and develop

The Tobii C12 is a full featured AAC device that gives users unlimited communication opportunities.

With the Tobii C12 users can communicate via text or symbols that generate synthesized or digitized speech, connect* with others through e-mail, text message or chat. They can write and edit text, surf the Internet, play music and games, control their environment, and manage their day – making independence easier than ever before.

Module based eye control

The Tobii C12 can be controlled in many ways including the built in touch screen, keyboard, mouse, switches, or the optional eye control unit, the Tobii CEye. Switching input methods is easy thanks to the modular construction.

The eye tracker unit docks seamlessly into the Tobii C12 and operates accurately regardless of glasses, contact lenses, eye color or light conditions. For those with uncontrollable head movements, the eye tracker's reliability

With its touch screen, weather resistant casing, hot swappable batteries, powerful speakers and optional eye control the Tobii C12 is the perfect choice for AAC users in any environment.

and large track box makes for effective eye control without the need for recalibration.

The communication software of choice

Users are free to choose the communication package that best fits their needs.

The Tobii C12 comes pre-installed with Tobii Communicator standard edition, which allows users to communicate using text or the more than 15,000 included SymbolStix® symbols. It comes with ready made page sets that help therapists set up an individual's language. The optional Communicator Premium supports e-mail, text messaging, and environmental control.

Easy to set up, use and personalise

The Tobii C12 is easy to setup and adjust to a user's input needs, voice preferences and other options. Multiple profiles can be stored, making it easy to select a profile according to an input change, language profile or to share a unit among users.

Personalise the Tobii C12 by replacing the side panels to give the device a discrete or colourful appearance. With the built-in camera it is possible to capture an image and turn it into a communication topic "on the fly".

Built to be used

The Tobii C12 is lightweight and can easily be carried. With many different mounting systems it is possible to mount the C12 on any wheelchair, walker or bed frame.

Thanks to its design, long lasting hot swappable batteries, powerful processor and stable operating system, this device offers users the strongest, fastest and most

flexible AAC device on the market. The Tobii C12 even allows for an active lifestyle thanks to its shock-mounted hard disk drive, silent operation and powerful stereo speakers.

The device has no buttons on the front and nothing protruding from the sides that can distract the user. The flat, moisture sealed surface makes for easy cleaning.

Product specifications Tobii C12

Display	12.1" resistive touch screen
Dimensions (WxHxD)	12.3" x 10.4" x 1.8" 312 x 264 x 46 mm
Weight	6,5 lbs 2,9 kg
Speakers	4 integrated
Microphones	2 integrated
Processor	IntelCore Duo U2500
RAM	2 GB DDR2
Operating System	Microsoft Windows 7
Hard disk drive	60 GB
SD Card	1 GB
Connectors	2 x USB 2.0, 1 x Micro USB, 1 x RJ45, 2 x switch, 1 x head phone (stereo), 1 x DC in
ECU	Integrated
Battery time (hot swappable)	5h
Camera	640x480 pixels
Deskstand	Integrated
Mounting systems supported	Optional brackets for Vesa, Daesyy and Rehadapt
Warranty	2 years

Accessories

The Tobii C12 comes bundled with all that is needed to get started. A list of accessories is available at www.tobii.com.

*Requires Tobii Communicator Premium and an internet connection.

© Tobii®. Illustrations and specifications do not necessarily apply to products and services offered in each local market. Technical specifications are subject to change without prior notice. All other trademarks are the property of their respective owners.

EUROPE, MIDDLE EAST, AFRICA

Tobii Technology AB
Karlsrovägen 2D
S-182 53 Danderyd
Sweden
Phone: +46 8 663 6990
Fax: +46 8 30 14 00
sales@tobii.com

CENTRAL EUROPE

Tobii Technology GmbH
Niedenau 45
D-60325 Frankfurt am Main
Germany
Phone: +49 69 24 75 03 40
Fax: +49 69 24 75 03 429
sales.de@tobii.com

NORWAY

Tobii Norge
Thormøhlens gate 55
Datablokken, 6 etasje
5008 BERGEN
Norway
Phone: +47 55 55 10 60
Fax: +47 55 55 10 61
sales.no@tobii.com

ASIA

Tobii Technology, Ltd.
3-4-13 Takanawa, Minato-ku
Tokyo 108-0074
Japan
Phone: +81-3-5793-3316
Fax: +81-3-5793-3317
sales.jp@tobii.com

NORTH AMERICA

Tobii ATl
333 Elm Street
Dedham, MA 02026
Toll-Free: 800-793-9227
Phone: 781-461-8200
Fax: 781-461-8213
Web: www.tobiati.com
sales.us@tobii.com

tobii
www.tobii.com